
BD125X/2


Rev. 19 / 06.2024

Bowie-Dick Test Pack

Result reference guide


Unprocessed

Unsatisfactory

Satisfactory

NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.

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ATTENTION: If your chemical indicator has this logo printed on it, it means that your product is compatible with exclusive **Trazanto** Automatic Quality Control and Traceability System for cleaning, disinfection and sterilization processes.
This System interprets the indicator's results in an accurate and automatic manner, and stores all the process data through Bionova® Software. Use **Trazanto** scanner to interpret the results (see instructions for use in the manual of the device).
ATTENTION: do not write on the indicator before scanning.

Indicators 100% Lead Free

Explanation of Symbols

Product designed for testing air removal in dynamic-air-removal steam sterilizers.

CE mark.

Authorized representative in the European Community.

Medical Device.

Batch number.

Expiration Date.

Manufacturer.

Bowie-Dick Test Pack

For vacuum-assisted steam sterilizers

Description

Chemdye® Bowie-Dick Test Pack detects air leaks, inadequate air removal and steam penetration level in vacuum-assisted steam sterilizers. The Chemdye® Bowie-Dick Test Pack consists of a chemical indicator sheet (125 mm x 125 mm) between porous material sheets and wrapped forming a package, with a steam indicator label.

The central sheet, with a regular pattern, corresponds to a Bowie-Dick classic indicator, i.e. BD Test Sheet. The chemical indicator changes from yellow to dark brown/black when processed. Any unexpected color change, such as lighter areas in the center of the sheet or different colors in the edges (i.e., a non-uniform color change) indicates an air pocket during the cycle due to sterilizer malfunction. Bowie-Dick Test Pack should be used daily, before running the first load of the day, after a sterilizer is installed or relocated, after a sterilizer malfunction, after sterilization process failures and after any major repairs of the sterilizer. A shortened cycle (i.e., a cycle omitting the post-vacuum drying phase) should be run first to properly heat up the sterilizer.

Instructions for use

1. Place the Bowie-Dick Test Pack horizontally with its label facing up in an empty vacuum-assisted sterilizer chamber, in the bottom/front of the sterilizer rack, and over the drain.
2. Run the cycle for no longer than 3.5 min. at 134 °C / 4 min. at 132 °C / 16.9 min. at 121 °C.

NOTE: If any serious incident occurs in relation to the device, it should be reported to Terragene S.A. and the competent authority of the State in which the user is established.

Interpretation of results

After completing the cycle, retrieve and examine the chemical indicator sheet.

Satisfactory test: It is indicated by a test sheet that shows completely uniform color change to dark brown/black. Such a test indicates rapid steam penetration, adequate air removal and lack of significant air leaks.

Unsatisfactory test: It is caused by an air pocket within the chamber and indicates incomplete air removal from the chamber, poor steam quality, superheated steam and/or non-condensable gases in the steam supply. If a non-uniform color change is noticed, it should be reported immediately to the supervisor and the sterilizer should not be used until examined for possible malfunction. Recall and quarantine the previous loads sterilized in that machine until it's certain that the machine has not failed before.

Warning

Do not use a Bowie-Dick Test Pack in gravity displacement steam sterilization cycles (where there is no pre-vacuum phase to check) or in other kinds of sterilization processes.

Storage

Store Chemdye® Bowie-Dick Test Packs in a dark place, 30-80 % RH at 10-30 °C.

It is recommended to protect the indicator Test Sheet after use inside a plastic foil for its storage.

Shelf life

Chemical Indicators have an expiration date of 5 years from the date of manufacture when stored at recommended storage conditions. Both, product expiration date and batch number are specified on product packaging and in each individual card. Do not use indicators after their expiration date.

Endpoint stability: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard along with cardboard waste according to your country's healthcare regulations.

The Bowie-Dick Test Pack is manufactured according to ISO 11140-1:2014 and ISO 11140-4:2007 standards.

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