
BT60 Biological Indicator

Rev. 12 / 11.2020

BIONOVA®

Spore Strips

Quality certification
Certificado de calidad
Bionova® BT60/

EO

STEAM

DRY

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Kill time / Tiempo de muerte

Z-value / Valor Z

Parameters determined at time of manufacture according to ISO 11138-1:2017, ISO 11138-2:2017, 11138-3:2017, 11138-4:2017 and IRAM 37102: 1999 (Parts 1, 2 and 3) standards. The values shown are reproducible only under the same conditions under which they were determined.

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ISO and USP Compliant
ATCC is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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80

30

EN Spore Strip
For Steam, Ethylene Oxide and Dry Heat sterilization

Composition
Each Bionova® BT60 envelope consists of a paper strip inoculated with a *Geobacillus stearothermophilus* ATCC 7953 and *Bacillus atrophaeus* ATCC 9372 spore population.

Product description
The envelope containing Bionova® BT60 Spore Strip is specifically

designed to control Steam (Steam), Ethylene Oxide (EO) and Dry Heat (DRY) sterilization processes. Spore strips should be incubated in an appropriate growth medium for the corresponding *Geobacillus stearothermophilus* ATCC 7953 or *Bacillus atrophaeus* ATCC 9372 growth. It is suggested the use of Bionova® MC20-2 Culture Medium for growth of *Geobacillus stearothermophilus* ATCC 7953 at 60 °C, MC1030-2 for growth of *Bacillus atrophaeus* ATCC 9372 at 37 °C and/or MC1020-2 for growth of both spores, selective according to the incubation temperature (*Geobacillus stearothermophilus* 60 °C / *Bacillus atrophaeus* 37 °C).

If the sterilization process was not successful, MC20-2, MC1030-2 or MC1020-2 Growth Medium will turn to yellow after incubation at the corresponding temperature, thus indicating the existence of living spores on the strips. If the sterilization process was successful, MC20-2, MC1030-2 or MC1020-2 Indicating Medium will remain its original color. A final readout should therefore be performed after 24 hours of growth medium incubation at 60±2 °C (MC20-2 and/or MC1020-2) for steam sterilization processes control, or 48 hours at 37±2 °C (MC1030-2 and/or MC1020-2) for Ethylene Oxide and Dry Heat sterilization processes control.

Precautions
Do not use Bionova® BT60 Spore Strip envelopes to control radiation cycles or other sterilization processes different from those for which they have been designed. Do not reuse spore strip envelopes.

Storage
Best stored in original package under the next conditions: temperatures between 10-30 °C, 30-80 % Relative Humidity.
Do not freeze.
Do not store near sterilizing agents or other chemical products.

Instructions for use
1. Identify the Bionova® BT60 envelope by writing the sterilizer number (in case of having more than one), load number and processing date on its back.
2. Pack it in an appropriate package along with material to be sterilized, according to recommended sterilization practices. Place the envelope in those areas you consider a priori more inaccessible to the sterilizing agent. A typical problematic area is the center of the load.
3. Sterilize as usual.
4. After the sterilization process has finished, you can:
A. Remove the Bionova® BT60 envelope from the package for processing and incubation.
B. For Ethylene Oxide sterilization, you can aerate the package containing the biological indicator with the rest of the load or remove the biological indicator from the sterilization package for incubation prior to material aeration.
5. Tear up the envelope at one end very carefully under sterility conditions (eg. laminar flow cabinet) and transfer spore strip with a sterile clamp to MC20-2, MC1030-2 or MC1020-2 Growth Medium tube, or to any other appropriate growth medium. **IMPORTANT:** Use latex gloves and chin strap when transferring spore strips from the envelope to growth medium tube. **WARNING!** Avoid spore strip contact with any kind of surface including the outer part of Growth Medium tube and of Bionova® BT60 envelope. **IMPORTANT:** Use a spore strip which has not been subjected to the sterilization process as a positive control each time a processed strip is incubated. Positive control guarantees that the incubation was carried out under appropriate conditions.
6. Incubate processed spore strips along with the strip used as positive control at the corresponding temperatures: 60±2 °C or 37±2 °C. It is advisable to make observations every 10 hours.
Color change of MC20-2, MC1030-2 or MC1020-2 Growth Medium to yellow means a failure in the sterilization process. If after incubation, no color change is visible on the growth medium, the result is negative (i.e. the sterilization process was successful). MC20-2, MC1030-2 or MC1020-2 Indicating Medium used for incubating the spore strips which have not been sterilized (positive control) should turn to yellow in order for the results to be valid. Record the positive results and discard them immediately as indicated below.
WARNING! Do not reuse the sterilizer until spore strip growth result is negative (MC20-2, MC1030-2 or MC1020-2 Growth Medium containing processed spore strip remains its original color).

Shelf life
Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date.

Disposal
Discard Growth Medium tubes already used to analyse Bionova® BT60 Spore Strips growth according to your country's healthcare and safety regulations. Positive Growth Medium tubes can be sterilized in gravity air displacement steam sterilizers at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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Terragene S.A