
BT97 Biological Indicator

Rev. 6 / 02.2024

Super Rapid Readout Fluorescence System

Quality certification Certificado de calidad Bionova® BT97 BioSurf

For airborne and surface disinfection by vaporized and aerosolized Hydrogen Peroxide / Para la desinfección de ambientes y superficies por vapor o aerosol de Peróxido de Hidrógeno. *Geobacillus stearothermophilus* ATCC® 7953

LOT

Population / Población _____ CFU / UFC

D - value / Valor D _____ sec./seg. (2 mg/l VH₂O₂, 50 °C)

Parameters determined at time of manufacture according to ISO 11138-1:2017 and IRAM 37102-1: 1999 standards. The values shown are reproducible only under the same conditions under which they were determined.

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ATCC® is a registered trademark of American Type Culture Collection.

Lic. Adrián J. Rovetto
Director Técnico
Technical Director

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Explanation of Symbols

VH202 Product designed for use with Vaporized Hydrogen Peroxide sterilization cycles
LOT Batch number

Manufacture Date
Expiration Date
Manufacturer

EN Biological Indicators

For airborne and surface disinfection by vaporized and aerosolized hydrogen peroxide

Composition

Bionova® BT97 BioSurf Biological Indicator consists of two tubes, A and B. Tube A contains a population of *Geobacillus stearothermophilus* ATCC® 7953 spores embedded in a stainless steel carrier, which is attached to the lid. Tube B contains a purple growth indicator medium within a glass ampoule and a chemical indicator of exposure in the external label.

Product description

Bionova® BT97 BioSurf Biological Indicator has been designed for quick and easy monitoring of airborne and surface disinfection processes by hydrogen peroxide.

The system consists of two tubes: tube A which contains the spores to challenge the disinfection process and tube B which is used to reveal its success.

Precautions

Do not use Bionova® BT97 BioSurf Biological Indicators to control processes different from the indicated in this instructions for use. Do not reuse the biological indicators.

Instructions for use

1. Inside the room to be disinfected open tube A and place the lid down and carrier up on the holder specially designed for the product. The rest of the tube may be discarded. Thus, *Geobacillus stearothermophilus* spores are exposed to evaluate the disinfection process. Place tube B without uncovering in the remaining position of the holder.

WARNING: Place the holder with the indicator in those areas that you consider a priori more inaccessible to the hydrogen peroxide. Places subject to minimal flows are suggested such as corners of the room, sectors in and around the equipment and spaces between disposable materials that will be used in the room.

To define the most difficult areas to decontaminate, it is recommended to use several indicators in the same process.

NOTE: A scheme of the room is attached to the product, where you can mark the position of the indicators with respect to the disinfection device.

2. Carry out the process of decontamination and aeration as usual.

3. Start the Bionova® Auto-reader Incubator, select the temperature (60 °C) and the required readout program (1 hour).

IMPORTANT: Always wear both sterile clothing and personal protective equipment (glasses, mask, gloves) to enter the enclosed area. Avoid sudden movements of air. Do not reveal Bionova® BT97 BioSurf Biological Indicator until the incubator reaches the required temperature.

4. Verify that the chemical indicator printed on the label of tube B turned to green. A color change confirm that the indicator has been exposed to hydrogen peroxide.

IMPORTANT: This color change is not evidence that the process was successful. If the chemical indicator color has not changed, review the process.

5. Crush the ampoule contained in the biological indicator with an individual ampoule crusher or with the ampoule crusher placed within incubator's incubation area.

Then shake the tube down vigorously, until the medium reaches the base of the tube and soaks the filter entirely. Incomplete wetting of the filter may lead to an incorrect fluorescence readout.

6. Remove the lid of tube B by taking the tube through the base. Discard it keeping the tube in vertical position.

7. Place the exposed lid of tube A on tube B so that the carrier contacts the culture medium.

8. Press the lid until the plastic tube is completely closed.

IMPORTANT: Use a non-processed biological indicator as a positive control every time a processed indicator is incubated. The positive control ensures that the incubation conditions were adequate.

9. Incubate the processed biological indicator along with the indicator used as a positive control for a maximum of 1 hour at 60 ± 2 °C for Super Rapid Readout.

Super Rapid Readout: 1 hour

Super Rapid Readout must be carried out in Bionova® auto-reader incubators. Final readout of negative results is readily available after 1 hour of incubation at 60±2 °C. Fluorescence readout is an indirect measure of the viability of *Geobacillus stearothermophilus* spores after the disinfection process (positive results). Furthermore, a failure in the process can also become evident by the color change of the medium. Due to the high sensitivity of the fluorescence results at 1 hour, conventional incubation for color change of Bionova® BT97 BioSurf Biological Indicator is not an advantage.

Visual confirmation: 48 hours

Optionally, you can perform a visual color-change confirmation after a 48-hour incubation. If the disinfection process has not been successful, culture medium will turn to yellow during incubation at 60 °C, indicating the presence of living spores. If disinfection was successful, culture medium will remain purple after the incubation process.

The positive control must show a color change from purple to yellow for results to be valid.

Visual readout: 7 days

A 7-day readout is optional and not intended to be routinely performed. This is an initial validation of the 1-hour readout. The 1-hour fluorescence results are compared to the 7-day visual readouts. The sensitivity of the system is determined as the difference between the 7 days positive indicators and false negative indicators (negative fluorescence readout and visually positive) in relation to the positives at 7 days. Based on the claimed ≥ 97 % readout reliability of the 1 hour biological indicator, there is no advantage in incubating the indicator beyond 1 hour.

NOTE: If 7-day readout is performed, a humidified environment will be required to prevent medium from drying out.

Storage

Store in a dark place and at temperatures between 10-30 °C, 30-80 % relative humidity. Do not freeze. Do not store biological indicators near sterilizing agents or other chemical products.

Shelf life

Biological Indicators have an expiration date of 2 years from the date of manufacture when stored at recommended conditions. Do not use indicators after their expiration date. Chemical Process Indicator on SCBI label have an expiration date of 2 years when used as part of SCBI. Endpoint Stability Reaction: chemical indicator endpoint shall remain unchanged for a period of not less than 6 months when stored at previously indicated conditions.

Disposal

Discard biological indicators after use according to your country's healthcare and safety regulations. The positive biological indicators can be autoclaved in a gravity air displacement steam sterilizer at 121 °C for 30 minutes, 132 °C for 15 minutes or 134 °C for 10 minutes; or in a dynamic air removal steam sterilizer at 132 °C for 4 minutes or 135 °C for 3 minutes.

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